TWISTER® Plus

Rotatable Retrieval Device

INSTRUCTIONS FOR USE

MODE D'EMPLOI

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INSTRUCTIONS FOR USE

DEVICE DESCRIPTION

The TWISTER[®] PLUS Rotatable Retrieval Device consists of a flexible wire/ catheter and a basket with net that is extended and retracted from the outer sheath using a three-ring handle. Device specifications are available as follows:

| Catheter Diameter | Catheter Length | Endoscope Minimum Working Channel |
|-------------------|-----------------|-----------------------------------|
| 2.5 mm | 230 cm | 2.8 mm |
| 1.7 mm | 120 cm | 2.0 mm |
| 1.7 mm | 160 cm | 2.0 mm |

INTENDED USE

The TWISTER[®] PLUS Rotatable Retrieval Device is intended to be passed through an endoscope's working channel to remove foreign body objects, such as polyps and food boluses.

WARNINGS

- Do not attempt to repair the device.
- For single use only. Do not reuse, reprocess or re-sterilize.
- Do not use if sterile barrier is damaged. If damaged, call Customer Service.
- Used product and packaging should be disposed of in accordance with hospital, administrative and/or local government policy.
- Instructions for Use should be read thoroughly as failure to follow the instructions for use may lead to serious clinical complications.
- Use of this device should be restricted to use by or under the supervision of physicians trained in endoscopic procedures.
- Reusing, reprocessing or re-sterilizing may compromise the structural integrity of the device and/or lead to device failure, which in turn may result in patient injury, illness or death. Re-use, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection including, but not limited to, transmission of infectious disease(s) from one patient to another.

PRECAUTIONS

- Inspect the TWISTER[®] PLUS Rotatable Retrieval Device and packaging for damage prior to use. Do not use product if opened or damaged. Confirm the product is consistent with the package label. Contact Customer to report and replace damaged product.
- The TWISTER[®] PLUS Rotatable Retrieval Device requires a working channel as referenced in the table specified in the Device Description section.
- Any use of the product, other than as indicated in the instructions, is not recommended.

PREPARING FOR USE

- Open and remove the TWISTER[®] PLUS Rotatable Retrieval Device from the pouch.
Note: Uncoil the TWISTER[®] PLUS Rotatable Retrieval Device before continuing with preparation to avoid damage.
- The TWISTER[®] PLUS Rotatable Retrieval Device should be visually inspected for kinks and complete product integrity. If any damage is observed, do not use the TWISTER[®] PLUS Rotatable Retrieval Device and contact Customer.
- The TWISTER[®] PLUS Rotatable Retrieval Device should be tested prior to being passed through the endoscope's working channel. Test the device by sliding the finger grip forward and backward. The flexible wire and basket with net should completely extend from and retract back into the tip of the catheter.
- Verify the working channel of the endoscope to ensure compatibility with the TWISTER[®] PLUS Rotatable Retrieval Device, (referenced in the table specified in the Device Description section).
- The TWISTER[®] PLUS Rotatable Retrieval Device is now ready for clinical use

INSTRUCTIONS FOR USE

Note: The TWISTER[®] PLUS Rotatable Retrieval Device should be advanced through the endoscope's working channel using short and precise movements to mitigate potential unintended damage to the device or endoscope.

- Identify the foreign body endoscopically and advance the TWISTER[®] PLUS Rotatable Retrieval Device through the endoscope's working channel until the distal end of the catheter is endoscopically visible.
- The distal end of the catheter should be advanced slightly distal to the foreign body however it is not recommended to advance the device if the foreign body has blocked the entire lumen. The TWISTER[®] PLUS Rotatable Retrieval Device should be opened completely by advancing the handle forward until it stops. Endoscopically confirm the net is fully open.
- The TWISTER[®] PLUS Rotatable Retrieval Device should be positioned over the foreign body and rotated into position by rotating the device handle. For best results, the foreign object should be centered or proximal in the net prior to closure.
- To secure the foreign body, retract the handle until the net closes around the foreign body. The net should be carefully opened and closed with gentle traction applied until the entire foreign body is retrieved. Excessive traction may cause damage to the device and/or the foreign body.
- The TWISTER[®] PLUS Rotatable Retrieval Device should then be retracted toward the endoscope so that it is visible, yet not compromising the endoscopic view to ensure visibility during removal of the device and foreign body.
- Once the scope, device and foreign body have been removed, the foreign body can be collected from the TWISTER[®] PLUS Rotatable Retrieval Device by advancing the handle forward to open the net. The foreign body should be processed per procedural standards.
- Once retrieval is complete, the device should be closed by advancing the handle and removed from the endoscope's working channel.

CONTRAINDICATIONS

- General medical condition that would not allow tolerance to endoscopy and/or other manipulation required.
- Patients demonstrating resistance.
- The device may not be used for applications other than for which it is intended.
- The device must not be passed through an incompatible endoscope.
- The device must not be re-sterilized and reused.
- The device is not intended to be implanted.

Any serious incident that has occurred in relation to the device must be notified to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Accessing Electronic Versions of IFU - Step-by-Step Guide

- Start by opening a web browser on your device and enter the official website URL for Horizons Corp: <https://www.horizonscorp.com>.
- Once the website loads, you will be directed to the main page, which features various tabs and options. Look for the tab labeled "Products."
- Click on the "Products" tab to access a dropdown menu that categorizes different product options. Within this tab, you can choose the specific product category that aligns with your requirements. Alternatively, you can enter the following URL: <https://www.horizonscorp.com/products> which will take you directly in the product pages
- Within your chosen product category, you will find a list of individual devices or models. Select the specific device that corresponds to your needs.
- After selecting the device, the main page will load with comprehensive information about the chosen product. At this point, you will have multiple options to explore, including "Specifications," "IFU," "Brochure," or "Video."
- To access the Instructions for Use (IFU), click on the "IFU" option. This will lead you to a section where the latest revision of the readable IFU is presented in an electronic format

7. In the IFU section, you will have the opportunity to view and read the most recent version of the IFU directly on the webpage. Additionally, if you prefer, you can choose to download a PDF copy of the IFU for offline reference.
8. If you are interested in accessing previous revisions of the IFU, you can find these versions within the same IFU section. Look for the category labeled "Previous Versions."

By following these steps, you will be able to access and utilize electronic versions of the Instructions for Use (IFU) for Horizons Corp products. Remember to make use of the IFU to ensure the proper and safe use of the products as intended.

| Symbol | Description | Symbol | Description |
|--------|---|--------|----------------------------------|
| | Country of manufacture | | Translation |
| | Date of manufacture | | Do not use if package is damaged |
| | Use by date | | Do not re-use |
| | Batch Code | | Consult instructions for use |
| | Catalogue number | | Important cautionary information |
| | Sterilized using ethylene oxide | | Do not re-sterilize |
| | Caution Federal (USA) law restricts this device to sale by or on the order of a physician | | Latex Free |
| | Keep Away from Sunlight | | Keep Dry |
| | Type BF applied part | | Medical Device |
| | Importer | | Authorized Representative |
| | Unique device identifier | | Sterile Barrier |

This product is protected by one or more of the following United States Patents: 9,101,342; 8,858,567.

TWISTER® is a registered trademark of Horizons International Corp.

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