



The **Power** of a Two in One
Rotatable **Snare** and Polyp **Catcher**

TORNADO™

One single instrument
One Single Procedure



Horizons International Corp.

We listen, Innovate, and Create to save lives


"The power of a Two in one"
Rotatable: Handle and tip
rotation is 1 to 1

TORNADO™

Rotatable, Snare and Polyp Catcher

- Rotatable, Snare and Polyp Catcher
- Unique, simple, advanced instrument
- Saves you time when you need it the most!
- Each instrument is packaged individually and supplied sterile
- For single use only

Step 1 SNARE




Ties enable full view of the polyp

Step 2 Cauterizer



Ties are released as the polyp is being cut.

Step 3 Semi Basket



Semi Basket is deployed to capture the polyp.

PRODUCT SPECIFICATIONS TORNADO™

ORDER NUMBER	CATHETER DIAMETER	TORNADO SIZE	LENGTH	BOX / OF	TYPE
TORNADO-O-23	2.5 mm	13 mm	230 cm	10/BOX	ROTATABLE / COLON
TORNADO-O-24	2.5 mm	22 mm	230 cm	10/BOX	ROTATABLE / COLON
TORNADO-O-25	2.5 mm	25 mm	230 cm	10/BOX	ROTATABLE / COLON

TORNADO™ is a Trademark of Horizons International Corp.

Protected by one or more of the following US Patents: 8,858,567 and 8,142,443.

Caution: Federal (USA) law restricts this device to be sold by, or on the order of a physician.

Distributed by:



We listen, Innovate, and Create to save lives

Horizons International Corp.

P.O. Box 213-3006 Zona Franca Metro,

Barreal de Heredia, Costa Rica

Phone: (506) 2261-9114 / Fax: (506) 2261-9115

FOR ORDERS IN THE U.S. PLEASE CALL TOLL FREE 866-245-7697

website: www.horizonscorp.com / **email:** info@horizonscorp.com