



ENDOSCOPE CLEANING BRUSHES

Horizons International Corp endoscope cleaning brushes are designed to facilitate disinfection and sterilization of scopes by removing organic matter and residue from internal channels and valve openings. Without proper cleaning, disinfection and sterilization of the scope cannot be achieved.

Choose from a complete selection of reusable and disposable cleaning brushes for gastrointestinal and pulmonary applications. Reusable brushes have controllable ring handles and flexible Teflon[®]-coated spring sheaths for easy clean-up. These brushes really stand up to hard wear and repeated use.

The 400 Series cleaning brushes are disposable and an economical alternative for endoscope cleaning. The 200 Series cleaning brushes are also disposable, and feature a spring-mounted brush for greater flexibility at the tip. A luer is attached to the proximal end for injection of the cleaning solution which flushes around the brush as it cleans. The brush has a blunt end cap to prevent possible damage to the scope channel and a built-in guidewire provides improved rigidity for insertion into the endoscope.

Horizons International Corp.

P.O Box 213-3006 Zona Franca Metro, Barreal de Heredia, Costa Rica
Phone (506) 2261-9114 • Fax (506) 2261-9115
www.horizonscorp.com • info@horizonscorp.com



ENDOSCOPE CLEANING BRUSHES

Reusable Scope Cleaning Brushes

Reusable scope cleaning brush with controllable ring handle and Teflon[®] coated spring sheath.

Product Number	Application	Brush Diameter	Recommended Minimum Channel Dia.	Working Length	Packaging Non-Sterile
Gastrointestinal					
	Pediatric Gastroscope	3mm	1.8 mm	160 cm	1 each
	Gastroscope/Colonoscope/Duodenoscope	4mm	2.8 mm	230 cm	1 each
	Gastroscope/Colonoscope/Duodenoscope	6mm	3.5 mm	230 cm	1 each
	Sigmoidoscope	6mm	3.2 mm	100 cm	1 each
Pulmonary					
	Bronchoscope	3mm	1.8 mm	115 cm	1 each

Disposable Scope Cleaning Brushes

The 200 series scope cleaning brushes feature a clear catheter and an injection port at the proximal end. The 400 series Blue Bullet[®] scope cleaning brushes feature a blue opaque tubing. Two specialty styles are available. The 478D (double end) has a cleaning brush attached to both ends of the tubing. The 478C (combination style) has a cleaning brush on one end and a scope valve cleaning brush on the other.

Product Number	Application	Brush Diameter	Recommended Minimum Channel Dia.	Working Length	Packaging Non-Sterile
Gastrointestinal					
	Gastroscope/Colonoscope	6 mm	3.5 mm	240 cm	50/box
	Pediatric Gastroscope	3 mm	1.8 mm	180 cm	50/box
	Gastroscope/Colonoscope/Duodenoscope	5 mm	2.8 mm	250 cm	50/box
	Pediatric Gastroscope	3 mm	2.0 mm	180 cm	50/box
	Combination valve/ channel cleaning brush for Gastroscope/Colonoscope/Duodenoscope	10 mm 5 mm	N/A 2.8 mm	230 cm	50/box
	Double-end cleaning brush for Gastroscope/Colonoscope/Duodenoscope	5 mm	2.8 mm	250 cm	50/box
Pulmonary					
	Bronchoscope	3 mm	1.8 mm	115 cm	50/box
	Bronchoscope	3 mm	1.8 mm	115 cm	50/box
Scope Valve Cleaning Brush					
	Endoscope/Bronchoscope	10 mm	N/A	5 cm	50/box

Horizons International Corp.

P.O Box 213-3006 Zona Franca Metro, Barreal de Heredia, Costa Rica
 Phone (506) 2261-9114 • Fax (506) 2261-9115
www.horizonscorp.com • info@horizonscorp.com