



360° Rotation
(1:1 Ratio In All Directions)

TWIST™ ROTATABLE SNARES

SINGLE-USE
ⓧ





Rotatable



Horizons International Corp.

We listen, Innovate, and Create to save lives

TWIST™ ROTATABLE SNARES



FEATURES	BENEFITS
Unique design	360° Rotation
Rotatable Snare	1:1 Ratio in all directions
Accuracy	Unique and simple to maneuver
Effectiveness	Saves you time during collection
Presentation	Each instrument is supplied in sterile packaging individually for single use only

PRODUCT SPECIFICATIONS TWIST™ ROTATABLE SNARES

PRODUCT CODE	CATETER DIAMETER MM	WORKING LENGTH CM	SNARE SIZE MM	TYPE	PACKAGING
HIC-SODRT-23	2.5	230	13	OVAL/HOT	BOX 10
HIC-SODRT-24	2.5	230	25	OVAL/HOT	BOX 10
HIC-SODRT-25	2.5	230	32	OVAL/HOT	BOX 10
HIC-SHDRC-23	2.5	230	10	HEXA/COLD	BOX 10
HIC-SODRC-23	2.5	230	10	OVAL/COLD	BOX 10

TWIST™ SNARE is Trademark of Horizons International Corp.

Caution: Federal (USA) law restricts this device to be sold by, or on the order of a physician.

Distributed by:



We listen, Innovate, and Create to save lives

Horizons International Corp.

P.O. Box 213-3006 Zona Franca Metro,
Barreal de Heredia, Costa Rica

Phone: (506) 2261-9114 / Fax: (506) 2261-9115

FOR ORDERS IN THE U.S. PLEASE CALL TOLL FREE 866-245-7697

website: www.horizonscorp.com / **email:** info@horizonscorp.com