


Rotatable Retrieval
Device

TWISTER® PLUS





Rotatable



Horizons International Corp.

We listen, Innovate, and Create to save lives

TWISTER® PLUS



FEATURES	BENEFITS
360° rotation (1:1 Ratio in all directions)	Designed to facilitate the capture and removal of polyps, food boluses and foreign bodies
Three Deep dimensional net	Engineered to capture and remove multiple polyps, large food boluses and foreign bodies in a single pass
Durability: interwoven net design	Net is attached to the outer wire using a interwoven design to help maintain the integrity of the net
Short throw handle	Ergonomically designed for nurse and technician comfort
Four net sizes	15 mm, 22 mm, 26mm and 32 mm net sizes are available to meet specific procedural requirements

PRODUCT SPECIFICATIONS TWISTER® PLUS

PRODUCT CODE	DESCRIPTION	Net Size (mm)	Sheath OD (mm)	Working Length (cm)	Minimum Working Channel (mm)	PACKAGING
HIC-GDL-27-5	TWISTER® PLUS Rotatable Retrieval Device	22	2.5	230	2.8	Box 5
HIC-GDL-27-10	TWISTER® PLUS Rotatable Retrieval Device	22	2.5	230	2.8	Box 10
HIC-GDL-26-5	TWISTER® PLUS Rotatable Retrieval Device	26	2.5	230	2.8	Box 5
HIC-GDL-26-10	TWISTER® PLUS Rotatable Retrieval Device	26	2.5	230	2.8	Box 10
HIC-GDL-29-5	TWISTER® PLUS Rotatable Retrieval Device	32	2.5	230	2.8	Box 5
HIC-GDL-29-10	TWISTER® PLUS Rotatable Retrieval Device	32	2.5	230	2.8	Box 10
HIC-GDL-30-5	TWISTER® PLUS Rotatable Retrieval Device	15	1.7	120	2.0	Box 5
HIC-GDL-30-10	TWISTER® PLUS Rotatable Retrieval Device	15	1.7	120	2.0	Box 10
HIC-GDL-31-5	TWISTER® PLUS Rotatable Retrieval Device	15	1.7	160	2.0	Box 5
HIC-GDL-31-10	TWISTER® PLUS Rotatable Retrieval Device	15	1.7	160	2.0	Box 10

TWISTER® PLUS is a registered Trademark of Horizons International Corp.
Protected by one or more of the following US Patents: 9,101,342; 8,858,567.

Caution: Federal (USA) law restricts this device to be sold by, or on the order of a physician.

Distributed by:



We listen, Innovate, and Create to save lives

Horizons International Corp.

P.O. Box 213-3006 Zona Franca Metro,
Barreal de Heredia, Costa Rica

Phone: (506) 2261-9114 / Fax: (506) 2261-9115

FOR ORDERS IN THE U.S. PLEASE CALL TOLL FREE 866-245-7697

website: www.horizonscorp.com / **email:** info@horizonscorp.com